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Instructions for use

INDICE / SUMMARY

ITA

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D.STEP CAST AND D.STEP CAST PLUS

EN

The cast walkers D.STEP CAST and D.STEP CAST PLUS are designed to be an **immobilizer for the foot and the ankle**, they allow to maintain constant the plantar contact surface and reduces the cutaneous stress.

The company FGP S.r.l. Unipersonale declares that the D.STEP CAST braces are medical devices designed and built according to the state - of - the art and they are made comply with the European Regulation 2017/745 in applicable parts.

It's important to consider, though, that **a misuse and/or any positioning by the patient without authorization and**

specific instructions given by the healthcare professional can lead to injuries or damages.

Before using the product, it is essential to read these instructions thoroughly.

Proper application is essential for the effective operation of the device.

During the D.STEP CAST and D.STEP CAST PLUS use, often check for redness and/or irritation on the skin. If you notice skin changes or feel discomfort contact the healthcare professional immediately.

The foot and leg when not bandaged, must be protected with a stocking. Be careful when walking on wet surfaces.

INDICATIONS

Brace designed to immobilize the foot and ankle.

Allows to maintain the surface of foot contact and limit skin stress.

Useful, for example, in cases of:

- Acute phase treatment and secondary prevention of diabetic foot ulcers (D.STEP CAST).
- Treatment of Charcot's foot (D.STEP CAST).
- Treatment of vascular foot.
- Treatment of fractures and amputations of the midfoot, forefoot and hindfoot (D.STEP CAST).
- Post operative treatment of the hindfoot lesions (D.STEP CAST PLUS).
- Acute phase treatment of the hindfoot ulcer (D.STEP CAST PLUS).

CONTRAINDICATIONS

No notes at the moment. The correct application according to the precise indications provided by the attending

healthcare professional, does not give rise to contraindications.

Application: for adults only.

LIMITATION OF LIABILITY

The company FGP S.r.l. Unipersonale assumes no responsibility for the use of the product in the conditions listed below as improper use:

- non-compliance with therapeutic indications;
- non-compliance with manufacturer's requirements and instructions;
- non-compliance or use contrary to current safety regulations;
- significant change in device size;
- significant alteration of the structure with modification, addition or elimination of components;
- product adaptation with significant shape alteration;
- sports, whether or not reduced;
- cleaning with unsuitable products;
- unsuitable storage;
- product abuse;
- driving vehicles and/or means of transport while using the product;
- assembly or disassembly of device components by unqualified personnel;
- failure to comply with assembly or disassembly instructions.

TECHNICAL CHARACTERISTICS

Diabetic foot orthosis made with new FGP Tecnocast technology.

STABILISING STRUCTURE

- Outer shell in high-rigidity technopolymer.
- Patented double-shell structure, Advanced Clamping Tech.
- Anti-slip rocker sole with boat-shaped design and provision for heel offloading wedges.
- Heel off-loading wedge, available in the height chosen by the healthcare professional (15mm or 30mm) (D.STEP CAST PLUS).

CLOSURE SYSTEM

- Equipped with 2 independent patented Spin Lock tension modulation systems (patented).

- Operable via rotation system for closure and pull-out for opening.
- Tensioners with ergonomic attachment and quick-release pull-tab for easy handling.

MODULAR PADDING SYSTEM

- Outer grippy material, inner open-cell foam for comfort and breathability.
- Additional EVA padding (30° Shore) for the calf, attachable via grippy system at the blue-petrol stitching area (optional, for volume compensation).

INTERNAL INSOLE

- Negative-pressure plantar system ensures total-contact with the plantar surface.
- Pneumatic circuit restores correct configuration each time the orthosis is worn.
- Insole is protected inside a pocket equipped with padded lateral protections for increased patient safety and comfort.

NOTES D.STEP CAST

Upon request, a set of heel off-loading wedges can be provided. The set consists of 2 wedges, 15 mm and 30 mm in height, which can be attached to the sole using dedicated screws in the area marked on the sole.

WARNINGS

- The product shall only be used on prescription or on the advice of an authorised health professional. Any adjustments and adaptations must be made only by healthcare professional (doctor or licensed health care professional).
- In order to ensure the optimal functioning of the pneumatic insole, each time the brace is removed for dressing changes, or at most every 3 days, the insole must be returned to its neutral position so that when it is repositioned, it adapts to the new position of the patient's foot. In the case of a patient authorized by the healthcare professional to remove the brace at home, it is recommended to reset the insole every day.
- Follow the cleaning and storage instructions carefully
- Be very careful when walking on slippery or wet ground.
- Avoid practices that make the product a tool or toy.
- Driving vehicles and/or means of transport' while using D.STEP CAST is not recommended.
- Do not expose the product to flames

MEASURE (AMBIDEXTROUS)

SIZE	SHOE SIZE
MEDIUM	37-42
LARGE	43-46

ORDER CODE:

DBT-DSTEPCAST

DBT-DSTEPCASTPLUS

- or sources of excessive heat.
- Make sure all closures are firm before use and walking.
- In the case of use in the heel lift configuration, make sure that this is well fixed to the sole before use.
- If wounds are present protect as per medical instructions.
- Use on intact skin.
- Do not use in case of wear, breakage or device's material integrity loss.
- The packaging and instructions are part of the medical device.
- Modification/alteration of the medical device excludes any liability of the manufacturer.
- To avoid a possible loss of functionality of the product or a limited duration of the product, it is recommended to prevent/limit contact between the device and creams, oils, lotions or fatty or acidic substances.
- Do not change the adjustments made by the healthcare professional on the medical device so as not to lead to a loss of functionality of the product.
- D.STEP CAST is designed to be used in combination with the DBT-RICCUNEI and DBT-DSTEPBALANCE devices made by FGP S.r.l. Unipersonale.
- D.STEP CAST PLUS can be used in combination with the DBT-DSTEPBALANCE device made by FGP S.r.l. Unipersonale.

INSTRUCTIONS FOR USE

CLEANING

RIGID FRAME

- Wash by hand using mild soap and rinse in cold water.
- Let it dry in the air or with a soft cloth.
- Avoid any aggressive liquid such as alcohol, bleach, ammonia, peroxide solvents and optical brighteners.
- Avoid using abrasive rubbers, brushes, steel wool, glass paper and knife blades for stain removal.
- Keep away from heat and sun.
- Do not use dryers or dryers.

PADDINGS

- Hand wash the paddings, the pads and the insole cover. It is possible to remove the insole cover from the cut located on the rear side.
NOTE: normal handwashing with water at max 30°C with neutral soap. Rinse thoroughly.
Let it dry in the air before use.
- Do not bleach.
- Do not dry in rotary drum tumble drier.
- Do not iron.
- Do not dry clean.

STORAGE

Store in a dry place at room temperature. Protect from heat, direct exposure to sunlight and moisture.

Disposal: do not disperse in the environment.

The product can be disposed of with normal household waste.

NOTE: This Class I medical device must be used with the appropriate warnings given in these instructions, in combination with those provided by the healthcare professional, which form an integral part thereof.

For any alleged problems arising as a result of the use of the device, contact the attending healthcare professional. In case of complaints or alleged malfunctions, please contact the retailer where you purchased the product. The orthosis is made for the exclusive use of a single patient. In the event of a serious accident, death or serious deterioration of the patient's health, related to the **MD**, contact the manufacturer and the competent authority immediately.

POSITIONING:


INSTRUCTIONS FOR THE HEALTHCARE PROFESSIONAL

In the case of positioning by the patient and/or their caregiver with the express authorization of their healthcare professional, follow the instructions and the training received from the healthcare professional in the previous months.

Unhook the padding elastics from the 4 fastening hooks.

Open each of the 6 tensioners and secure it to its base.

Pull the knobs outward to release tension on the pods.

After opening the sides, remove the padding.

Put on the padding, making sure that your heel is positioned at the back. Close the straps so that the padding is properly fitted. Verify that there is no overlap of material, wrinkles or pressure points.

Close first the right "wing" and then the left one, paying attention that the front surfaces of the latter overcome those of the right side, behind the teeth of the left side must be placed in the respective housings.

With the brace open, position the foot so that the back of the leg meets the mast of the hull.

Close the belts and after pushing in the knobs put the pods in slight tension without tightening.

Insert the tube supplied with the electric pump into the connector located at the rear of the hull.

At this point tighten the tensioners (pod) creating the desired closure around the leg of the patient.

With the electric pump at hand, ask the patient to stand up, placing their weight on the leg wearing the walker, and then take a few steps.

At this point, activate the electric pump by pressing the green button. The air suction will stop automatically, then disconnect the tube from the pneumatic circuit.

Let the patient walk for a few steps and verify that the brace does not move, if necessary tighten the tensioners further. Check that the patient's heel does not touch the edge of the brace to avoid the risk of hyper pressure. If the patient has no muscle tone in the calf, use the pads that can be positioned using the hook-and-loop fastening system supplied with the device.

Attach the rubber bands present on the padding on the appropriate housing of the hull.

To return the insole to its neutral position, press and hold the pneumatic circuit valve button located at the rear of the hull to allow air to re-enter the system. At the same time, use your other hand to evenly redistribute the balls inside the sole, restoring its uniform configuration.

At the end of the positioning phase, the healthcare professional is required to communicate to the patient the two values shown on the graduated scales, determined by the overlap of the two sides, so that, in the event of repositioning requested and authorized by the healthcare professional, the patient is able to correctly replicate the tensioning of the PODs.

WARNING: in order to ensure the optimal functioning of the pneumatic insole, each time the brace is removed for dressing changes, or at most every 3 days, the insole must be returned to its neutral position so that when it is repositioned, it adapts to the new position of the patient's foot.

In the case of a patient authorized by a healthcare professional to remove the brace at home, it is recommended to replace the sole every day.

NOTES FOR THE HEALTHCARE PROFESSIONAL

If the patient has no muscle tone in the calf, use the pads that can be positioned using the hook-and-loop fastening system supplied with the device.

The following items are included:

- a rear pad for volumetric filling, supplied pre-applied to the padding; remove it only if deemed unnecessary;
- two tibial pads and two instep pads, which can be used either for filling or for protection if the patient experiences pain in these specific areas specifiche zone.

The rear pad is attached directly to the padding, following the outline indicated by the stitching for correct positioning.

If necessary, soft pads must be positioned by the healthcare professional and must not be tampered with by patients. For this reason, they are attached directly to the wings of the brace. To position them, remove the fabric rectangles protecting the hook strips, then apply the pads to create the appropriate pressure relief areas. In addition, please take special care to ensure that the pads do not interfere with the closure of the brace. The tibial and instep pads can be cut to size as needed.

TECHNICAL INSTRUCTION

HEEL WEDGES MOUNTING

In case an off-loading of the patient's heel is necessary, D.STEP CAST is equipped with two wedges that can be attached to the sole that are able to move the patient's load towards the forefoot.

1 Remove the four pre-cut sole parts.

2 This will expose four threaded inserts where the wedge for the heel off-load will be placed.

3 Choose the wedge of the appropriate height, based on the therapeutic needs, place it so that the highest part is located at the heel.

- 4 Fix the wedge to the hull with the dedicated screws. Each wedge is equipped with two pairs of screws of different lengths, the longer screws serve to fix the heel area, the shorter ones the area towards the forefoot. No special tightening torque is required.

- 5 Verify that the wedge is correctly positioned so that there is no distance between the surfaces in contact and that there is no movement between the parts.

- 6 In D.STEP CAST PLUS whenever is necessary to change the wedge, proceed bottom up as described to take off the old wedge and, then, place the new one.

PNEUMATIC SYSTEM REPLACEMENT

In the event of a malfunction of the pneumosphere system that can be related to a puncture of the insole or a malfunction of the non-return valve, it is suggested to replace the entire system (cod. DBT-RICPIS).

- 1** Unscrew the two screws holding the pneumatic circuit protection cover in place.

- 2** Remove the cover.

- 3** Cut the tube that arrives to the pneumatic circuit.

- 4** Remove the tube from its housing by removing the insole from the hull and replacing it with the new one.

- 5** When replacing the insole, check that the tube is correctly positioned in its housing in the hull.

- 6** Connect the tube to the new pneumatic circuit.

- 7** Reposition the pneumatic circuit on the appropriate housings.

- 8** Reposition the protective cover and secure it with the appropriate screws and washers. Test the circuit.

INSTRUCTIONS FOR THE PATIENT

WARNING: the removal and repositioning of the device by the patient must only take place after an appropriate written request, authorization, and instruction provided by the healthcare professional.

During each check-up, the professional is required to inform the patient of the tension to be applied to the PODs, indicating the value on the graduated scales that corresponds to the overlap of the sides.

REMOVAL OF THE DEVICE

The device must only be removed by the patient upon written request, authorization, and adequate instruction provided by the healthcare professional. The instructions for removal are as follows:

1. Unhook the padding elastic bands from the 4 fastening hooks.
2. Pull the knobs outward to release tension on the pods.
3. Open each of the 6 tensioners and secure it to its base.
4. Open the sides and remove the brace.

POSITIONING OF THE DEVICE

The device must only be positioned by the patient upon written request, authorization, and adequate instruction provided by the healthcare professional.

For repositioning, refer to the section 'Instructions for the healthcare professional'.

WARNING: Before repositioning the device, the insole must be returned to its neutral position so that it adapts to the new position of the foot.

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Consultare le
istruzioni d'uso
Refer to instructions
for use

Conservare in
ambiente secco
Store in dry
environment

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AZIENDA CERTIFICATA
CERTIFIED COMPANY
ISO 13485:2016

Scarica l'App
FGPOrteSystem®

FISTRDSTEFPCAST
Rev.01 - 02/26