D.STEP
Total Contact Diabetic Cast

A TOTAL CONTACT INNOVATION

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A technological
revolution
grounded
in science.

FGP
Orthopedic
& Rehabilitation
Equipment

FGP E MEDI	5
<i>FGP AND MEDI</i>	
TOTAL CONTACT CASTING	6
<i>TOTAL CONTACT CASTING</i>	
IL FUTURO. D.STEP, TOTAL CONTACT DIABETIC CAST	12
<i>THE FUTURE. D.STEP, TOTAL CONTACT DIABETIC CAST</i>	
TEST DI LABORATORIO	19
<i>LABORATORY TESTS</i>	
SCHEDA TECNICA	22
<i>TECHNICAL DATA SHEET</i>	
RICAMBI	24
<i>SPARE PARTS</i>	
ACCESSORI	25
<i>ACCESSORIES</i>	
BIBLIOGRAFIA	26
<i>REFERENCES</i>	

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FGP and MEDI. Two protagonists of well-being with a common vision.

FGP is a company that brings Italian manufacturing excellence and a passion for personal well-being to the world. We develop orthoses, rehabilitation devices, and orthopaedic solutions with an integral approach to quality. Our commitment is Made in Italy, built on a network of qualified partners and professionals, with the aim of providing everyone with the freedom to move forward again in safety and comfort.

We believe in a future where everyone can face mobility challenges with confidence and peace of mind. We strive to become the preferred partner of physicians, orthopaedic technicians, therapists, and healthcare professionals, turning technology and artisan know-how into an effective, comfortable, and lasting recovery experience.

At the beginning of 2023, we became part of the international MEDI group, one of the world's leading manufacturers of compression stockings and medical products. A union born from the sharing of values and a common vision, which has consolidated our know-how and strengthened our position among the most innovative companies in the production of orthopaedic supports or orthoses. Together, we continue to carry forward Italian excellence with even more strength and international credibility.

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Cfr. [6], [19], [20], [22].

Overview of Diabetic Foot Ulcers

Neuropathic ulcers are skin lesions that result from a combination of factors. Some of the main factors that can lead to DFUs include:

- 1 Peripheral neuropathy** (damage to the peripheral nervous system);
- 2 Shear forces** (typically friction between the skin and footwear);
- 3 Repetitive stress** (resulting, for example, from walking);
- 4 Pathological peak pressure**, associated with anatomical alterations that create areas of overload while standing;

These factors rarely cause problems in healthy individuals, but sensory neuropathy secondary to diabetes leads to the loss of pain sensitivity. As a result, the natural protective mechanism – pain-induced off-loading of the affected limb – is diminished or even lost. The result is that foot injuries in diabetic neuropathic patients are often underestimated or ignored due to the absence of pain. The patient continues to walk, despite the injury, thereby perpetuating the pathological process that caused the ulcer to appear.

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Cfr. [8], [19], [22].

Off-loading techniques

Obtaining an effective off-loading system is not straightforward, particularly because the period of time required for the wound to heal is often lengthy.

Many patients need to maintain a certain level of mobility throughout the day, as prolonged bed rest is not a practical option for extended periods. Other off-loading techniques – such as wheelchairs or crutches – have not received validation from proven scientific data.

See. [2], [15].

THE OFF-LOADING METHODS CURRENTLY IN USE CAN BE FUNCTIONALLY CLASSIFIED INTO TWO MAIN GROUPS:

REMOVABLE DEVICES, which allow the patient to control how long they are worn;

NON-REMOVABLE DEVICES, which are not removable by the patient and must be removed or changed by healthcare personnel.

Several reviews on off-loading devices have shown that non-removable total-contact casts are the most effective option for off-loading DFUs.

Ankle joint immobilisation reduces pressure in the forefoot and plantar regions during the propulsive phase of gait. Both total-contact casts and rigid removable walkers can immobilise the ankle. However, casts ensure consistent use, regardless of patient compliance.

The total-contact cast shortens stride length, slows foot motion, and decreases the forces applied to the foot. It has been measured that the total-contact cast is able to reduce plantar pressure by 32%, 63% and 69% on the fifth, fourth, and first metatarsal heads, respectively; by 65% on the big toe; and by 45% on the heel.

A total-contact cast can be used to off-load the rearfoot but there is no agreement among experts that this approach is as effective as the one for releasing pressure on the forefoot.

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Cfr. [3], [4], [12].

The **TOTAL CONTACT CAST (TCC)** is a specially designed plaster cast developed to reduce peak plantar pressures by redistributing loads across a wider surface area (off-loading) in patients with ulcerative lesions (DFUs) resulting from neuro-ischaemic complications of diabetes mellitus. The reduction of pressure on the lesion has been shown to be an essential therapeutic component in the healing of plantar ulcers affecting the midfoot and rearfoot.

The use of a Total-Contact Cast (TCC) requires specific expertise, manual skill, and time to prepare. The cast, usually made of fiberglass, encloses the foot (including the toes) and extends up the leg to the mid-calf. This configuration enables the redistribution of plantar loads onto the leg, thereby relieving pressure at the ulcer site.

This reduction of loads removes the traumatic component from the lesion allowing a more physiological and rapid process of tissue regeneration. The TCC eliminates flexion-extension of the ankle during walking, thereby reducing the shear and torsional forces that constitute an additional factor that hinder the tissue repair of the lesion.

Effective offloading is the key element in the management of diabetic foot syndrome ulcers, particularly in neuropathic plantar lesions of the forefoot and midfoot. Sensorimotor neuropathy contributes to ulcer formation by both altering foot anatomy –creating areas of overload at the metatarsal heads, midtarsal region, toe tips, and heel – and disturbing gait biomechanics, leading to prolonged weight-bearing on the metatarsophalangeal regions.

Together with the assessment and management of infections and vascular status, TCC is an essential and irreplaceable component in the treatment of diabetic ulcerative lesions, being the most effective and reliable method for off-loading DFUs.

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Clinical evidence

The use of the total-contact cast is the method that clinical evidence indicates as the most effective for the relief of pathological pressure. Randomised clinical trials and meta-analyses have provided strong evidence on the increase in the healing rate and reduction in the re-epithelialisation times of DFUs.

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Armstrong et al. (2001) demonstrated that TCC achieved a healing rate of 89.5%, compared with 65% for a removable walker and 58.3% for footwear. The average healing time was 4.8 weeks with the TCC, 7.2 weeks with the removable walker, and 8.7 weeks with therapeutic footwear. These findings were confirmed by Lavery more than 14 years later in a new study.

Lavery et al. (2014) demonstrated that TCC achieves a healing rate of 88.9%, compared with 50% for acute-phase orthoses, and 40% for orthoses designed to control shear stress. The mean healing time was 5.4 weeks with TCC, 8.9 weeks for acute-phase orthoses, and 6.7 weeks with shear-stress-control orthoses.

Figure 1


Figure 1 - Randomized, controlled clinical trial from 2014, comparing total contact casting (TCC), healing sandal (HS) and shear-reducing walking boot (SRB). Data replotted from Lavery et al. 2014.[15]

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Contraindications

TCC is contraindicated in specific clinical conditions, as reported in the related scientific literature. The presence of untreated infection, osteomyelitis or in patients with critical ischaemia who have not undergone revascularisation, TCC is not recommended by the related literature (these subjects are even excluded from clinical trials with TCC). Following the state of the art, it should be used with caution in cases of deep or heavily exuding wounds, and in ataxic, blind, or severely obese patients.

Patient compliance is necessary in minimising complications. TCC ensures compliance, hence it's a non-removable device. However, use of the device may reduce the patient's quality of life, as they may require assistance with daily activities, such as using the toilet or getting in and out of bed. These limitations lower TCC's acceptability. In addition, the use of advanced dressing techniques or therapies requiring daily clinical monitoring may not be compatible with the use of a non-removable device. In addition, TCC carries a moderate risk of causing new friction-related lesions and is not recommended for claustrophobic patients.

Prolonged use of TCC has been documented to cause joint stiffness and muscle atrophy, due to the biomechanical role of the device. It can also worsen significant postural instability in susceptible patients. It's acknowledged, TCC effectiveness strongly depends on its correct application. Highly specialized technicians is needed to ensure the correct construction of the device. This personnel is often not available in the reference clinical spots (it's reported only the 2-3.7%). Moreover, the technician and the constitutive material are a sensitive cost, since TCC must be periodically built (once per week for minimum 6-8 weeks) and patients show high recidivism (it's estimated 20-58% of the subjects).

See [1], [2], [6], [15], [16], [23].

Charcot neuro-osteoarthropathy

Charcot Neuro-Osteoarthropathy (CN), as a complication of diabetic neuropathy, is a pathology that involves the progressive alteration of the anatomy of the foot and/or ankle, resulting in severe deformities with a high risk of ulceration, infection and, ultimately, amputation.

The pathophysiological process begins with an inflammatory phase and progresses to structural deformity. Areas of increased pressure contribute to the development of ulcerative lesions, which carry a high risk of amputation if an effective medical and/or surgical treatment program is not initiated. In the literature, the prevalence of ulcerative lesions is described as being 40-50% in patients with CN.

TCC is considered the gold standard not only for the treatment of DFUs, but also to protect the foot and ankle in the early stages of CN when the risk of bone and joint fractures and dislocations is high. It is widely accepted that the CCT should be regarded as the reference method for immobilisation and protection of acute-phase CN, as it is adapted, custom-fitted, and non-removable.

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An innovative structure for total-contact immobilisation.

D.STEP IS A BRACE FOR TIBIOTALAR (ANKLE) JOINT IMMOBILISATION CAPABLE OF CONTROLLING PLANTAR LOAD AND PROVIDING TOTAL CONTACT WITH THE PATIENT'S LIMB.

FGP has revolutionised the design, conceiving the splint as a single piece, with no removable parts, equipped with the innovative A.C.T. closure system, and a pneumatic system (P.I.S.) capable of replicating the plantar contour, thereby ensuring immobilisation, off-loading, and reduction of skin friction.

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TECNOCAST MAKES D.STEP A REVOLUTIONARY SPLINT WITH TOTAL-CONTACT TECHNOLOGY. IT FEATURES CUTTING-EDGE DESIGN SOLUTIONS (INCLUDING THE ADVANCED CLAMPING TECH PATENT) COMBINED WITH THE USE OF TECHNOLOGICAL MATERIALS TO PROVIDE A SOLUTION THAT IS TRULY INNOVATIVE AMONG ORTHOSES FOR THE DIABETIC FOOT.

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PATENTED CLOSURE-TENSION MODULATION SYSTEM THAT ALLOWS PRECISE ADJUSTMENT OWING TO:

- a rotational drive mechanism that easily allows the correct tension to be obtained;
- a tension-release mechanism activated by a simple “pull-out” movement of the knob;
- two differentiated tensioning circuits, for the upper part and for the lower part of the splint, allowing adaptation according to the patient’s specific needs.

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PATENTED “CLAMP” CLOSURE SYSTEM, EXCLUSIVE TO D.STEP, WHICH PROVIDES BIOMECHANICAL LOCKING OF THE TIBIOTALAR JOINT THANKS TO:

- patented double-shell structure;
- high-rigidity technopolymers;
- shell geometry designed to adapt to any anatomical conformation.

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PLANTAR SYSTEM WITH NEGATIVE-PRESSURE INSOLE, WHICH ENSURES TOTAL CONTACT WITH THE PLANTAR SURFACE THANKS TO:

- pneumatic circuit that utilises the patient’s weight to shape the insole for total-contact with the foot;
- system with vacuum pump that maintains the foot’s contour for prolonged wear;
- ability to reset and remodel the insole whenever necessary;
- lateral containment structure of the insole, which further protects the patient’s foot from any shear stress and medial and lateral marginal impingements.

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1 HIGH-RIGIDITY TECHNOPOLYMER DESIGN

Thanks to its characteristics and design, it stabilises and locks of the tibiotalar joint and foot, reducing friction.

2 WIDE MALLEOLAR OFF-LOADING

with contoured edge surfaces that prevent localized pressure ulcers.

3 LARGE POSTERIOR HEEL OFF-LOADING

with contoured edges to prevent localized pressure ulcers.

4 SPIN LOCK SYSTEM

allows precise adjustment of the brace closure tension.

5 MODULAR PADDING SYSTEM

adapts the brace to the anatomical conformations of the patient and stabilises the foot to prevent internal movement. Thanks to the attachment system on the structure, movements between the two elements are also prevented.

6 ANTI-SLIP ROCKER OUTSOLE

The boat-shaped design, featuring curves and containment edges, is engineered to provide correct off-loading, stability, and smooth movement.

7 LONGITUDINAL OFF-LOADING

designed to provide the necessary adaptability and rigidity for the patient.

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VACUUM INSOLE SYSTEM release button.

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System with **VACUUM PUMP** that maintains the foot's contour for prolonged wear.

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*DESPITE THE RIGOUR OF ITS
FUNCTIONALITY, ATTENTION WAS PAID TO
THE DESIGN OF THE PRODUCT SO THAT IT
WAS ALSO AESTHETICALLY PLEASING.*

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*TENSIONERS WITH EASY ATTACHMENT
SYSTEM facilitated by a pull-tab and
ergonomic shapes for easier handling.*

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*MODULAR PADDING system to ensure
maximum adaptability even in situations
of significant deformity. The main
component is made of internal open-cell
foam to ensure breathability and comfort.*

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Laboratory tests have demonstrated that the D.STEP brace is comparable to a total-contact cast in terms of plantar pressure distribution, both in static and dynamic conditions.

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

FOR THE COMPARISON, A D.STEP BRACE (FIGURE A) AND A TOTAL-CONTACT CAST APPLIED BY SPECIALIST PERSONNEL (FIGURE B) WERE USED.



ORTOSTASI / ORTHOSTASIS

GAMBALETTO GESSATO
TOTAL CONTACT CAST

Pressione
Pressure



TUTORE D.STEP
D.STEP BRACE

Area di contatto cm²
Contact area cm²



DEAMBULAZIONE / GAIT

GAMBALETTO GESSATO
TOTAL CONTACT CAST

Pressione
Pressure

TUTORE D.STEP
D.STEP BRACE

Area di contatto cm²
Contact area cm²

Misure (ambidestro)
Measures (bilateral)

TAGLIE SIZE	CALZATA SCARPA SHOE FIT
MEDIUM	37-42
LARGE	43-46

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Indications

Brace for immobilisation of the foot and ankle.

Designed to control plantar contact pressure and reduce cutaneous stress. Useful, for example, in cases of:

- Acute-phase treatment and secondary prevention of diabetic foot ulcers
- Charcot Foot Treatment
- Vascular foot treatment
- Treatment of fractures and amputations of the midfoot, forefoot and hindfoot

Technical features

DIABETIC FOOT BRACE MADE WITH NEW FGP TECNOCAST TECHNOLOGY

STABILISING STRUCTURE

- Outer shell in high-rigidity technopolymer
- Patented double-shell structure, Advanced Clamping Tech
- Anti-slip rocker sole with boat-shaped design and provision for heel offloading wedges

CLOSURE

- Equipped with 2 independent patented Spin Lock tension modulation systems (patented)
- Operable via rotation system for closure and pull-out for opening
- Tensioners with ergonomic attachment and quick-release pull-tab for easy handling

MODULAR PADDING SYSTEM

- Outer grippy material, inner open-cell foam for comfort and breathability
- Additional EVA padding (30° Shore) for the calf, attachable via grippy system at the blue-petrol stitching area (optional, for volume compensation)

INTERNAL INSOLE

- Negative-pressure plantar system ensures total-contact with the plantar surface
- Pneumatic circuit restores correct configuration each time the brace is worn
- Insole is protected inside a pocket equipped with padded lateral protections for increased patient safety and comfort

NOTES

Upon request, a set of heel off-loading wedges can be provided. The set consists of 2 wedges, 15 mm and 30 mm in height, which can be attached to the sole using dedicated screws in the area marked on the sole.

D.STEP CAST PLUS

TOTAL CONTACT DIABETIC CAST PLUS

CODICE D'ORDINE / ORDER CODE: DBT-DSTEPCASTPLUS

Misure (ambidestro)
Measures (bilateral)

TAGLIE SIZE	CALZATA SCARPA SHOE FIT
MEDIUM	37-42
LARGE	43-46

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Indications

Brace for immobilisation of the foot and ankle. Designed to control plantar contact pressure and reduce cutaneous stress. Useful, for example, in cases of:

- Acute phase treatment of hindfoot ulcers
- Vascular foot treatment
- Postoperative treatment of hindfoot lesions

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Technical characteristics

DIABETIC FOOT BRACE MADE WITH NEW FGP TECNOCAST TECHNOLOGY

STABILISING FRAME

- Outer shell in high-rigidity technopolymer
- Patented double-shell structure, ADVANCED CLAMPING TECH
- Anti-slip rocker sole with boat-shaped design and provision for heel offloading wedges
- Heel off-loading wedge, available in the height chosen by the clinician (15mm or 30mm)

CLOSURE SYSTEM

- Equipped with 2 independent patented Spin Lock tension modulation systems (patented)
- Operable via rotation system for closure and pull-out for opening
- Tensioners with ergonomic attachment and quick-release pull-tab for easy handling

MODULAR PADDING SYSTEM

- Outer grippy material, inner open-cell foam for comfort and breathability
- Additional EVA padding (30° Shore) for the calf, attachable via grippy system at the blue-petrol stitching area (optional, for volume compensation)

INTERNAL INSOLE

- Negative-pressure plantar system ensures total-contact with the plantar surface
- Pneumatic circuit restores correct configuration each time the brace is worn
- Insole is protected inside a pocket equipped with padded lateral protections for increased patient safety and comfort

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Padding made from open-cell foam to ensure breathability and comfort.

Complementary padding is included to guarantee optimum fit, even in cases of significant deformity.

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Extra thickness made from EVA foam, Shore 30°, attachable to the padding via the grippy system to provide volumetric compensation for the calf.

MISURE / MEASURES

M - L

CODICE D'ORDINE / ORDER CODE

DBT-RICPAD

SISTEMA COMPLETO VACUUM INSOLE SYSTEM COMPLETE VACUUM INSOLE SYSTEM

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Complete replacement of the pneumatic system consisting of:

- Vacuum-moulded insole with cover
- Pneumatic circuit with cap and screws

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Spare part consisting of:

- Vacuum pump with connecting tube
- USB charging cable

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Set of two wedges, 15 mm and 30 mm in height, made from technopolymer and attachable to the sole using dedicated screws. Equipped with anti-slip rubber at the base. These wedges can be optionally used with the brace to relieve pressure on the heel.

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Ambidextrous elevation system made of rubber, with non-slip sole. Equipped with a set of 3 EVA shims useful for correcting the height discrepancy between the two lower limbs and allowing for balanced walking.

MISURE MEASURES	M	L	XL
CALZATA SCARPA SHOE FIT	38-43	43-46	46-50
LUNGHEZZA LENGTH	28,5 cm	30,6 cm	34 cm

CODICE D'ORDINE / ORDER CODE

DBT-STEPBALANCE

Bibliografia / References

1. Bastiaan P. Vierhout, R. V.-P. (2022). Comparing a non-removable total contact cast with a non-removable softcast in diabetic foot ulcers: A retrospective study of a prospective database. *Diabetes Research and Clinical Practice*.
2. Carlo Caravaggi, E. F. (2000). Effectiveness and Safety of a Nonremovable Fiberglass Off-Bearing Cast Versus a Therapeutic Shoe in the Treatment of Neuropathic Foot Ulcers: A randomized study. *Diabetes Care*.
3. Carlotta Gauna, F. R. (2024). Offloading systems for the treatment of neuropathic foot ulcers in patients with diabetes mellitus: a meta-analysis of randomized controlled trials for the development of the Italian guidelines for the treatment of diabetic foot syndrome. Panel of the Italian Guidelines for the Treatment of Diabetic Foot Syndrome.
4. Caroline E. Fife, M. J. (2014). Diabetic Foot Ulcer Off-loading: The Gap Between Evidence and Practice. Data from the US Wound Registry. *Advances in Skin & Wound Care*.
5. Claire M. Capobianco, J. J. (2010). The role of an extended medial column arthrodesis for Charcot midfoot neuroarthropathy. *Diabetic Foot & Ankle*.
6. David G. Armstrong, H. C. (2001). Off-Loading the Diabetic Foot Wound: A Randomized Clinical Trial. *Diabetes Care*.
7. David G. Armstrong, L. A. (2004). It's Not What You Put On, but What You Take Off: Techniques for Debriding and Off-Loading the Diabetic Foot Wound. *Clinical Infectious Diseases*.
8. Fischer, D. H. (2022). Ventral Windowed Total Contact Casts Safely Offload Diabetic Feet and Allow Access to the Foot. *Journal of Diabetes Science and Technology*.
9. Greenhagen, R. a. (2009). Total contact casting for neuropathic ulcers: a lost art? *J Diabetic Foot Complications. The Journal of Diabetic Foot Complications*.
10. Isabella Gigante, E. D. (2023). Offloading of diabetes-related neuropathic foot ulcers at Swedish prosthetic and orthotic clinics. *Diabetes/Metabolism Research and Reviews*.
11. Jacqueline J. Wertsch, M., Lawrence W. Frank, M., Hongsheng Zhu, P., & Melvin B. Price, D. P. (1995). Plantar pressures with total contact casting. *Journal of Rehabilitation Research and Development*.
12. K. Jones, M. B. (2023). Rehabilitation for people wearing offloading devices for diabetes-related foot ulcers: a systematic review and meta-analysis. *Journal of Foot and Ankle Research*.
13. K. Morona, J. (2013). Comparison of the clinical effectiveness of different off-loading devices for the treatment of neuropathic foot ulcers in patients with diabetes: a systematic review and meta-analysis. *Wiley Online Library*.
14. Kurtis L Kowalski, J. D. (2002). Evaluation of Fiberglass versus Plaster of Paris for Immobilization of Fractures of the Arm and Leg. *Clinical Trial*.
15. Lawrence A. Lavery, K. R. (2014). Randomised clinical trial to compare total contact casts, healing sandals and a shear-reducing removable boot to heal diabetic foot ulcers. *International Wound Journal*.
16. Lawrence A. Lavery, S. A. (1996). Reducing Dynamic Foot Pressures in High-Risk Diabetic Subjects With Foot Ulcerations: A comparison of treatments. *Diabetes Care*.
17. Lewis J, L. A. (2013). Pressure-relieving interventions for treating diabetic foot ulcers (Review). *Cochrane Database of Systematic Reviews*.
18. Michael J. Mueller, D. R. (2003). Effect of Achilles tendon lengthening on neuropathic plantar ulcers: A randomized clinical trial. *The Journal of Bone and Joint Surgery*.
19. Piaggese, A. (2007). An Off-the-Shelf Instant Contact Casting Device for the Management of Diabetic Foot Ulcers: A randomized prospective trial versus traditional fiberglass cast. *American Diabetes Association*.
20. Raspovic, A. (2014). A survey of offloading practices for diabetes-related plantar neuropathic foot ulcers. *Journal of foot and ankle research*.
21. Rogers, L. C. (2011). The Charcot Foot in Diabetes. *Diabetes Care*.
22. Snyder, R. J. (2014). The management of Diabetic Foot Ulcers Through Optimal Off-Loading. *Journal of the American Podiatric Medical Association*.
23. Velasco-Rodríguez-Rabadán, S. (2023). Effectiveness of Off-Loading Devices in Patients With Active Diabetic Foot Ulcer: A Systematic Review. *The International Journal of Lower Extremity Wounds*.
24. Volker E. Schmidt, J. H. (1973). Mechanical Properties of Orthopedic Plaster Bandages. *Journal of Biomechanics*.
25. Wee Teck Lim, H. R. (2020). The Real-Life Outcome of VACOPed Boot in the Management of Diabetic Foot Ulcers. *The International Journal of Lower Extremity Wounds*.

AZIENDA CERTIFICATA
ISO 13485:2016

Scarica l'App
FGPOrteSystem®

fgpsrl.it

FGP srl Via Alessandro Volta 3, I-37062 Dossobuono VR
T +39 0458600867 F +39 0458600835 E fgp@fgpsrl.it

FGP

Orthopedic
& Rehabilitation
Equipment